


# The right Pipette Tip



Pipette tips Filter tips Volume range	Transferpette® single channel Nominal volume											Transferpette® multichannel Nominal volume														
	1 µl	2.5 µl	5 µl	10 µl	10 µl	20 µl	20 µl	25 µl	50 µl	100 µl	200 µl	200 µl	250 µl	500 µl	1000 µl	2 ml	5 ml	10 ml	10 µl	20 µl	25 µl	50 µl	100 µl	200 µl	300 µl	
0.1 - 20 µl	✓	✓ <sub>2)</sub>		✓ <sub>2)</sub>															✓	✓						
0.5 - 20 µl		✓		✓		✓													✓	✓						
1 - 50 µl		✓ <sub>2)</sub>		✓		✓													✓	✓						
2 - 200 µl <sup>1)</sup>			✓		✓		✓	✓	✓	✓	✓										✓	✓	✓	✓	✓	✓
5 - 300 µl			✓		✓		✓	✓	✓	✓	✓										✓	✓	✓	✓	✓	✓
50 - 1000 µl <sup>1)</sup>												✓	✓	✓	✓											
0.5 - 5 ml																✓	✓									
1 - 10 ml																		✓								
0.1 - 1 µl	✓	✓ <sub>2)</sub>		✓ <sub>2)</sub>															✓							
0.5 - 10 µl		✓		✓		✓													✓	✓						
1 - 20 µl				✓		✓													✓	✓						
2 - 20 µl			✓		✓		✓	✓	✓	✓	✓	✓									✓	✓	✓	✓	✓	✓
5 - 100 µl			✓		✓		✓	✓	✓	✓	✓	✓									✓	✓	✓	✓	✓	✓
5 - 200 µl											✓										✓	✓	✓	✓	✓	✓
50 - 1000 µl													✓	✓												

✓ = Tip volume less than pipette's nominal volume

<sup>1)</sup> Bulk tips are yellow or blue colored, racked tips are colorless in a yellow-colored or blue-colored tip tray

<sup>2)</sup> Calibration and possible adjustment needed